

Continuous Performance Test (CPT) for ADHD

Patient Information

Name:

Age:

Gender:

Date of Birth:

Date of Test:

Referring Physician:

Instructions

Please carefully read the instructions before starting the test. This test aims to assess your attention and impulsivity. You will be presented with a series of visual and auditory stimuli. Your task is to respond promptly to the target stimuli while ignoring non-target stimuli.

Test of Variables of Attention (TOVA)

Visual Test:

- The patient will be presented with geometric shapes on the screen.
- They are instructed to click the designated button or switch when they see the target shape.
- They should refrain from clicking when they see a non-target shape.

Auditory Test:

- The patient will hear sounds through headphones or speakers.
- They should press the designated button or switch when they hear the target sound.
- They should not respond when they hear a non-target sound.

Test Duration:

- Visual Test: Approximately 10.9 minutes for children ages 4-5, and 21.6 minutes for ages 6 and up.
- Auditory Test: Same duration as the visual test.

Scoring:

- Correct Responses: Number of accurate responses to target stimuli.
- Omission Errors: Number of missed responses to target stimuli.
- Commission Errors: Number of incorrect responses to non-target stimuli.
- Reaction Time: Average time taken to respond to target stimuli.

Results

Correct Responses:

Omission Errors:

Commission Errors:

Reaction Time:

Interpretation

Correct Responses:

Omission Errors:

Commission Errors:

Reaction Time:

Clinical Recommendations**Note to Healthcare Practitioner**